

Drybath90

Dry Incubator

Technical specifications

Setting temp. range	0°C~49.9°C(Ambient temp. must be lower than setting temp.)
Accuracy (control temp.)	≤±0.5°C
Accuracy (display temp.)	≤±0.1°C
Heating-up time	≤20min(Ambient temperature 10°C~30°C)
Shortcut key of setting temp	25°C,30°C,37°C
Input voltage	19V 2.1A
External dimension	185mm*175mm*85mm(length*width*height)
Weight	2Kg
Dimension	185*175*85mm(L*W*H)
Incubation Positions	22 positions total 15 positions for test tubes or direct reading cuvettes 5 positions for coagulation cups 2 positions for reagent bottles



Drybath90

Dry Incubator

Multifunctional Incubator, Smart Design



Shenzhen East Medical Technology Co.,Ltd
Add:190# longping Road,Huate Industrial Zone,longgang District,Shenzhen ,GuangDong China 518116
Tel : 0755-33275355 : Fax:0755-22658200
www.east-medical.com